

# M186

Renewal

## Sakamoto Whole Body Pregnant Model

● Accessories/External genitalia (skin color, transparent), Fetus (Delivery), Lubricant oil, Baby powder

### Useful for serial training of birth assist!



#### Specifications

Height	154cm
Weight	About 9kg

#### Exchangeable parts


- Exchangeable genitalia (skin color)
- Exchangeable genitalia (transparent)
- Exchangeable umbilical cord (10 pieces/set)

#### Practice


### Free-style delivery possible



▲Supine position



▲Squat position



▲Standing position



▲Lateral position



▲All fours (Knee-chest position)

Emergency +

Nursing


Injection

Experience


Acupuncture

Feature

1 Can see the status of fetal head through the transparent genitalia.



▲Protection of the perineum



▲Crowning



▲Delivery of posterior shoulder

2 Can check placenta/amnion



3 Can practice umbilical cord ligation and cutting



▲Exchangeable umbilical cord

M186-5

Sakamoto Whole Body Pregnant Model (with Heart sound Fetus)

●Accessories/External genitalia (skin color), Heart sound Fetus



M186-6

Fetus (Heart sound)

●Size/W26×D12×H15cm

●The actual heartbeat sound comes.



M186-3

Fetus (Delivery) (with placenta)

●Size/W40×D20×H10cm



M186-7

Leather Fetus (with placenta)



M186-8

Fetus (Palpation)

●Accessories/Exclusive base

